

# SARS-CoV-2 Neutralizing antibody Test Kit (Colloidal Gold)



## Basic Parameters

Specimen: whole blood, serum or plasma.  
 Test time: 10 mins  
 Packing: 1 Test/Box  
 Storage: 10-30°C (50-86°F)  
 Shelf life: 1 year.


**Accurate, efficient and reliable.**

## Purpose

SARS-CoV-2 Neutralizing antibody is an important marker for assessing the effectiveness of the new coronavirus vaccines. This reagent is suitable for neutralizing antibody detection in samples from individuals with vaccine injection people and recovered people.

## Significance


SARS-CoV-2 Neutralizing antibody Test Kit (Colloidal Gold) is of great significance for the development of new coronavirus vaccines, evaluation of effectiveness, and evaluation of neutralizing antibody levels in the population.



## Advantage of Product

- CE Marked
- Whole blood
- High specificity no false negative
- High sensitivity  $\geq 95\%$
- Rapid test without analyzer

## Principle of Test



## Interpretation of Result

